

SHOWERBUDDY SB2T TUBBUDDY TILT





SKU: SB2T

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The SB2t TubBuddy Tilt bath transfer system is designed to allow a caregiver to transfer the user safely and securely over the bathtub, offering a significantly cheaper alternative to a full bathroom remodel.

The tilt function has the additional benefit of providing pressure relief as well as reducing caregiver burden by allowing the user to transfer into the bathtub without having to raise their legs.

The SB2t can also roll over most standard toilets, or function as a standard toilet commode.

FEATURES

- · Rolls over most standard toilets
- Arms removable / lockable / flip-back
- Height adjustable
- No rust aluminium frame (Grade 6061 T6)
- · Removable and swing away footrests
- Safety belts lap & chest belts plus calf straps
- Caster Wheels 125mm x 4 locking
- Tool-free assembly

SPECIFICATIONS

- Overall Height 108cm-118cm
- Overall Depth 144cm (when footrests extended and fully tilted)
- Overall Width 57cm



- Seat Height 63cm-73cm
- Seat Depth 43cm
- Seat Width 47.5cm
- Tilt Angle 0-30 degrees (infinite positioning)
- Toilet Clearance 40cm-50
- Weight Capacity 250kg

WILL THIS FIT YOUR BATHROOM?

<u>Use Showerbuddy's Bathcheck app to see which Showerbuddy chair will work for you - Click Here!</u>

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 319903 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au

