



SHOWERBUDDY CUSHION CLOSED FRONT - SB6C/SB6W



SKU: SCCF-H

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

When seeking additional comfort and privacy during showering and commode use, the Showerbuddy Cushion Closed Front (SCCF) presents as a perfect add-on accessory for the Roll-InBuddy series users. This cushion is meticulously designed to seamlessly blend with the Roll-InBuddy SB6C/W models, ensuring an enhanced user experience.

FEATURES

- **Closed Front Design:** This feature offers users an extra layer of privacy by providing a contained commode gap, ensuring they feel secure and comfortable during use.
- **Tailored Fit:** Exclusively crafted for the Roll-InBuddy SB6C/W models, the SCCF cushion ensures a snug fit, eliminating concerns of mismatch or instability.
- **Durable Construction:** Made from PU foam, the same reliable material used in other Showerbuddy cushions, the SCCF promises longevity and resilience, even with regular usage.
- **Hassle-Free Cleaning:** Keeping the SCCF clean is a breeze. Its smooth surface ensures that it can be wiped down quickly, maintaining hygiene with minimal effort.
- **Versatile Installation:** Whether you've already invested in a Roll-InBuddy unit or are considering purchasing one, the SCCF can be easily added at any point, ensuring adaptability to users' changing preferences.

SPECIFICATIONS

- **Length:** 50.8 cm
- **Width:** 29.21 cm
- **Height:** 8.94 cm
- **Weight:** 0.9 kg
- **Material:** Closed-cell foam

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au