

SHOWERBUDDY OPEN FRONT CUSHION FOR SB6 HIGH PROFILE



SKU: SCSF-6H

PRODUCT DESCRIPTION

Purpose & Benefits

The **ShowerBuddy Open Front Cushion for SB6 High Profile** is a high-profile cushion accessory designed for use with the Roll-InBuddy SB6c and SB6w shower chairs. It features 50mm thick high density foam to provide enhanced cushioning and comfort. The cushion has a hand-stretched PU fabric covering with no edge join and is designed for easy cleaning. It is screw-fixed to the chair and comes with commode rails factory-fitted.

Features

- 50mm thick high density foam construction.
- Hand-stretched PU fabric covering with no edge join.
- Screw-fixed installation to the chair.
- Designed for use with the Roll-InBuddy SB6c and SB6w.
- Easy-to-clean PU fabric covering.
- Commode rails are factory-fitted.

Specifications

Specification	Details
Product Code	SCSF-6 H
Foam & Fabric	Soft - PU Skin
Foam Thickness	5 cm
Dimension A	18 cm
Dimension S	12 cm
Dimension T	9.5 cm
Overlay Cushion	SCOC
Insert	No

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
 1300 633 468

support@medimarnt.com.au

