




### SHOWERBUDDY PETITEBUDDY

**\$819.00 - \$889.00** Price range: \$819.00 through \$889.00



**SKU:** P1-H

## PRODUCT DESCRIPTION

### Purpose & Benefits of the ShowerBuddy PetiteBuddy

The ShowerBuddy PetiteBuddy is an accessory designed to adapt standard ShowerBuddy transfer systems for smaller users, including children and petite adults. By modifying the chair's dimensions, the PetiteBuddy ensures a comfortable and secure fit, enhancing safety and usability during bathroom transfers. This adaptability allows the ShowerBuddy system to grow with the user, providing a cost-effective solution that eliminates the need for frequent equipment upgrades.

### Features

- **Adjustable Backrest and Neck Pillow:** Offers customisable support to accommodate the user's specific comfort and positioning needs.
- **Reduced Width Between Armrests:** Narrows the space between armrests to better fit smaller body frames, ensuring a secure and comfortable seating position.
- **Smaller Commode Opening:** Features a reduced commode opening size, providing added safety and comfort for petite users during transfers.
- **Tool-Free Assembly:** Designed for easy installation without the need for tools, allowing for quick setup and adjustments as needed.
- **Compatibility:** The PetiteBuddy accessory is compatible with ShowerBuddy models SB1, SB2, SB2T, and SB3T, offering versatility across various transfer systems.

### Specifications

- **Backrest Width:** 380 mm
- **Backrest Depth Adjustment:** 325 mm to 375 mm
- **Width Between Armrests:** 420 mm

- **Commode Opening Size:** 170 mm
- **Compatibility:** Compatible with ShowerBuddy models SB1, SB2, SB2T, and SB3T
- **Assembly:** Tool-free assembly for easy installation