




## SHOWERBUDDY ROLL-IN BUDDY XXL SB6C-26

**\$1,387.00**



**SKU:** SB6c-26

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Showerbuddy Roll-In Buddy XXL SB6C-26 is a specialised shower chair designed to provide safe and convenient bathroom solutions for larger individuals or those with limited mobility. This model is part of the Showerbuddy range of bathroom mobility aids, known for their quality, safety, and ease of use. The Roll-In Buddy XXL is specifically engineered to accommodate the needs of bariatric users, ensuring comfort and security during bathroom activities.

### FEATURES

- **Bariatric Design:** The chair is built to support a higher weight capacity, making it ideal for bariatric users. Its sturdy frame ensures safety and stability.
- **Roll-In Capability:** Designed to easily roll into a shower, the SB6C-26 eliminates the need for transfers, reducing the risk of falls and injuries associated with bathroom activities.
- **Comfortable Seating:** The chair features a wide, comfortable seat with padding, providing adequate support and comfort during use.
- **Adjustable Components:** Comes with adjustable footrests and armrests, allowing for a customizable fit to meet the specific needs of each user.
- **Durable Construction:** Made with high-quality, corrosion-resistant materials suitable for use in wet environments.
- **Safety Features:** Equipped with locking casters to secure the chair in place during transfers and use, and a safety belt for added security.
- **Ease of Cleaning:** The materials used are easy to clean and maintain, ensuring hygiene.

## **SPECIFICATIONS**

- Chair Width – over caster wheels 670mm
  - Chair Width – over arms 760mm
  - Chair Width – inside arms 660mm
  - Chair Length – including footrest 970mm (730mm with footrests removed)
  - Backrest Width 650mm
  - Backrest Height from Cushion to top 380mm
  - Toilet Height clearance (maximum) 475mm
  - Toilet Width clearance at narrowest 560mm (at casters)
  - Toilet clearance can depend upon the shape of the toilet pan
  - Seat Cushion Dimensions 650mm W x 430mm D
  - Seat Height from Floor 510mm to 585mm (Adjustable)
  - Seat Cushion Opening 175mm W x 335mm L
  - Box Dimensions L 590mm x W 660mm x H 490mm
  - Box Gross Weight 17.5kg
- 

## **WILL THIS FIT YOUR BATHROOM?**

[Use Showerbuddy's Bathcheck app to see which Showerbuddy chair will work for you - Click Here!](#)