



## SHOWERBUDDY ROLL-INBUDDY SOLO

**\$1,103.00**



**SKU: SB6w**

## PRODUCT DESCRIPTION

### Purpose & Benefits

The **Roll-InBuddy Solo** allows users to regain their independence and comfort in the shower. With two 24" drive wheels, this self-propelled shower commode chair enables users to easily manoeuvre in and out of roll-in showers. Constructed from lightweight, durable aluminium, it supports up to 136kg and includes an 18" wide easy-clean cushion seat for comfort. The integrated commode system makes it a practical choice for home or facility use.

### Features

- 136kg Safe Working Load (SWL)
- Rolls over most standard toilets for added convenience
- Durable, low-maintenance, lightweight aluminium frame
- Lockable, flip-back armrests for easy transfers
- Height adjustable frame for customised positioning
- Removable, swing-away, and height-adjustable footrests
- No-rust aluminium with stainless steel parts and hardware

- Includes safety lap or chest belts plus a calf strap for security
- 100mm locking castor wheels for stability
- Commode system with deflector splash guard
- Soft blow-mould cushion for user comfort
- Tool-free assembly for quick setup

## Specifications

Specification	Detail
Chair Width Over Wheelchair Wheels	635 mm
Chair Width Over Arms	565 mm
Chair Width Inside Arms	470 mm
Chair Length	1060 mm (with footrest) / 820 mm (without footrest)
Backrest Width	448 mm
Backrest Height	380 mm (cushion to top)
Max Toilet Height Clearance	495 mm
Seat Cushion	460 x 430 mm
Seat Height From Floor	570 - 595 mm
Seat Cushion Opening	175 x 335 mm