



SHOWERBUDDY ROLL-INBUDDY SOLO



SKU: SB6w

PRODUCT DESCRIPTION

Purpose & Benefits

The **Roll-InBuddy Solo** allows users to regain their independence and comfort in the shower. With two 24" drive wheels, this self-propelled shower commode chair enables users to easily manoeuvre in and out of roll-in showers. Constructed from lightweight, durable aluminium, it supports up to 136kg and includes an 18" wide easy-clean cushion seat for comfort. The integrated commode system makes it a practical choice for home or facility use.

Features

- 136kg Safe Working Load (SWL)
- Rolls over most standard toilets for added convenience
- Durable, low-maintenance, lightweight aluminium frame
- Lockable, flip-back armrests for easy transfers
- Height adjustable frame for customised positioning
- Removable, swing-away, and height-adjustable footrests
- No-rust aluminium with stainless steel parts and hardware

- Includes safety lap or chest belts plus a calf strap for security
- 100mm locking castor wheels for stability
- Commode system with deflector splash guard
- Soft blow-mould cushion for user comfort
- Tool-free assembly for quick setup

Specifications

Specification	Detail
Chair Width Over Wheelchair Wheels	635 mm
Chair Width Over Arms	565 mm
Chair Width Inside Arms	470 mm
Chair Length	1060 mm (with footrest) / 820 mm (without footrest)
Backrest Width	448 mm
Backrest Height	380 mm (cushion to top)
Max Toilet Height Clearance	495 mm
Seat Cushion	460 x 430 mm
Seat Height From Floor	570 - 595 mm
Seat Cushion Opening	175 x 335 mm

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,

Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC**Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au