




## **SHOWERBUDDY SB1 TRANSFER SHOWER COMMODE**

**\$3,026.00**



**SKU: SB1**

## **PRODUCT DESCRIPTION**

### **PURPOSE & BENEFITS**

ShowerBuddy SB1 can help you gain access to your home shower again to enjoy a private shower in comfort and with dignity by allowing you to effortlessly glide over that shower step or threshold, seated comfortably for the entire transfer.

The ShowerBuddy SB1 shower transfer system will turn your existing bathroom into an accessible bathroom that allows you to enjoy showering again without the excessive cost, delays, and stress of a full bathroom remodel. With only one transfer, you can once again enjoy a shower in your own bathroom. The ShowerBuddy SB1 alleviates heavy lifting for family members and caregivers by eliminating multiple transfers during toileting and showering.

### **FEATURES**

- Adjustable height: The seat height can be adjusted from 550mm to 675mm to accommodate users of different heights.
- Removable arms: The arms can be removed for easier transfers or for storage.
- Rotates 360 degrees: The shower base can be rotated 360 degrees to allow the user to face the shower controls or to allow an assistant to help with cleaning.
- Safety belts: The chair is equipped with safety belts to prevent the user from falling.
- Footrest: The footrest provides support for the user's feet.
- Tool-free assembly: The Showerbuddy SB1 can be assembled without the need for tools.
- Lockable caster wheels

## **SPECIFICATIONS**

- Chair Width – over caster wheels 520mm
- Chair Width – over arms 570mm
- Chair Width – inside arms 475mm – smaller users see accessory pages
- Chair Length – including footrest 800mm (640mm excluding)
- Backrest Width - 460mm
- Backrest Height from Cushion to top - 440mm
- Toilet Height Clearance (maximum) - 510mm
- Toilet width clearance at narrowest - 295mm (centre of track, approx. 200mm above floor)
- Toilet clearance can depend upon the shape of the toilet pan
- Bridge Length (Standard) 405mm or Accessory 200mm or 300mm or customized, see Accessories
- Seat Cushion Dimensions 460mm W x 430mm D
- Seat Height from Floor 550mm to 675mm (Adjustable)
- Seat Cushion Opening 175mm W x 335mm L
- Shower Base unit dimensions 5000mm x 455mm
- Box dimensions 760mm x 570mm x 630mm
- Box Gross Weight 29.7kg

---

## **WILL THIS FIT YOUR BATHROOM?**

[Use Showerbuddy's Bathcheck app to see which Showerbuddy chair will work for you - Click Here!](#)