

SHOWERBUDDY SB3T IN-BUDDY TILT



SKU: SB3T

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Showerbuddy SB3T rolling shower chair is perfect for an accessible or roll-in wet area showers. The SB3T rolls over most toilets and offers a tilt that releases pressure for the user and improves comfort. The chair has a standard front opening cushion, head rest, swing away and removable footrest, flip back and removable armrest, 4" lockable casters, removable commode, seat belts, and height adjustable.

FEATURES

- Weight capacity – 160kg.
- Tilt angle 0-30 degrees, infinite positioning
- Rolls over most standard toilets
- Low maintenance
- Arms – removable / lockable / flip-back
- Height adjustable
- Footrests – removable and swing-away, height-adjustable
- No rust aluminium frame (Grade 6061 T6)
- Stainless steel parts and hardware (Grade 304)
- Safety belts – lap & chest belts plus calf strap
- Caster Wheels – 125mm x 4 locking
- Commode system
- Neck rest, fully adjustable
- Tool-free assembly

SPECIFICATIONS

- Chair Width – over caster wheels 490mm
- Chair Width – over arms 570mm
- Chair Width – inside arms 475mm – smaller users see accessory pages
- Chair Length – including footrest 1060mm (940mm with footrests fitted inwards)
- Safe Working Load - 160kg
- Chair Tilt 0-30 degrees with infinite positioning and lock
- Backrest Width 460mm
- Backrest Height from Cushion to top 440mm
- Toilet Height clearance (maximum) 535mm
- Toilet Width clearance at narrowest 375mm (at casters)
- Toilet clearance can depend upon the shape of the toilet pan
- Seat Cushion Dimensions 460mm W x 430mm D
- Seat Height from Floor 510mm to 610mm (Adjustable)
- Seat Cushion Opening 175mm W x 335mm L
- Box Dimensions L 840mm x W 580mm x H 760mm
- Box Gross Weight 27.4kg

WILL THIS FIT YOUR BATHROOM?

[Use Showerbuddy's Bathcheck app to see which Showerbuddy chair will work for you - Click Here!](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au