




## SHOWERBUDDY SB6C ROLL-INBUDDY LITE SHOWER COMMODE



**SKU:** SB6c

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Showerbuddy SB6c Roll-InBuddy Lite is a versatile and easy-to-use bathing and toileting aid for people with mobility impairments. It is designed to fit over most standard toilets and roll-in showers, making it a convenient and practical solution for a variety of bathroom needs.

The Roll-InBuddy Lite features a lightweight, rust-resistant aluminium frame that is easy to move and manoeuvre. It also has a comfortable padded seat and backrest, as well as adjustable footrests and armrests for added support. The seat height is adjustable, so it can be customised to fit the individual user.

For safety, the Roll-InBuddy Lite comes with four locking caster wheels and two safety belts. The caster wheels can be locked in place to prevent the unit from moving during use. The safety belts can be worn around the user's waist and chest to provide additional support.

The Roll-InBuddy Lite is also equipped with a commode system that can be used for toileting. The commode system is easy to install and remove, and it includes a deflector splash guard to prevent spills.

The Showerbuddy SB6c Roll-InBuddy Lite is a valuable bathing and toileting aid for people with mobility impairments. It is easy to use, safe, and versatile, making it a convenient and practical solution for a variety of bathroom needs.

### FEATURES

- **Weight capacity of 136kg:** The Roll-InBuddy Lite can support users up to 136kg, making it a

great option for people of all sizes.

- **Rolls over most standard toilets:** The Roll-InBuddy Lite has four locking caster wheels that make it easy to move around the bathroom. It can also roll over most standard toilets, so you can easily transfer from the toilet to the shower.
- **Lockable / flip-back arms:** The arms of the Roll-InBuddy Lite can be locked in place or flipped back, depending on your needs. This makes it a versatile option that can be used for both bathing and toileting.
- **Height adjustable:** The seat of the Roll-InBuddy Lite is adjustable, so you can find the perfect height for your comfort. This is especially important for people with limited mobility, as it allows them to get in and out of the shower or toilet more easily.
- **Removable and swing-away, height-adjustable footrests:** The footrests of the Roll-InBuddy Lite can be removed or swung away, depending on your needs. This makes it easy to get in and out of the shower or toilet, and it also provides additional space for you to move around.
- **No rust aluminium frame (Grade 6061 T6):** The frame of the Roll-InBuddy Lite is made from rust-resistant aluminium, making it a durable and long-lasting product.
- **Stainless steel parts and hardware (Grade 304):** The parts and hardware of the Roll-InBuddy Lite are made from stainless steel, which is rust-resistant and easy to clean.
- **Safety belts - lap or chest belts plus calf strap:** The Roll-InBuddy Lite comes with two safety belts: a lap belt and a chest belt. These belts can be worn to provide additional security while bathing or toileting.
- **Caster Wheels - 100mm all locking:** The four caster wheels of the Roll-InBuddy Lite are all locking, which provides added stability and safety.
- **Commode system with deflector splash guard:** The Roll-InBuddy Lite comes with a commode system that can be used for toileting. The commode system is easy to install and remove, and it includes a deflector splash guard to prevent spills.
- **Tool-free assembly:** The Roll-InBuddy Lite is easy to assemble, so you can get started using it right away.

## SPECIFICATIONS

- Chair Width – over caster wheels 505mm
- Chair Width – over arms 565mm
- Chair Width – inside arms 470mm
- Chair Length – including footrest 890mm (630mm with footrests removed)
- Backrest Width 448mm
- Backrest Height from Cushion to top 380mm
- Toilet Height clearance (maximum) 490mm
- Toilet Width clearance at narrowest 365mm (at casters)
- Toilet clearance can depend upon the shape of the toilet pan
- Seat Cushion Dimensions 460mm W x 430mm D
- Seat Height from Floor 510mm to 585mm (Adjustable)

- Seat Cushion Opening 175mm W x 335mm L
- Box Dimensions L 580mm x W 570mm x H 500mm
- Box Gross Weight 15.8kg

## **SIZING BROCHURE**

[Roll-InBuddy Lite SB6c - Sizing Brochure](#)

---

## **WILL THIS FIT YOUR BATHROOM?**

[Use Showerbuddy's Bathcheck app to see which Showerbuddy chair will work for you - Click Here!](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## **SA**

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## **VIC**

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## **NT**

**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)