

## SHOWERBUDDY SB7E ECO TRAVELLER



**SKU:** SB7e

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The ShowerBuddy SB7e Eco Traveller folding shower commode is lightweight and designed to be perfect for use in tight spaces, or for taking with you whilst travelling / storing away when not in use.

Fully height adjustable, the SB7e can roll straight over most toilets and into the shower, greatly reducing manual handling for carers, and simplifying transfers. The SB7e can also be used as a commode chair, with the bucket being removed and attached from the rear.

### FEATURES

- Weight capacity – 120kg
- Rolls over most standard toilets
- Arms – fixed height
- Adjustable Footrests – removable and swing-away, height-adjustable
- No rust aluminium frame (Grade 6061 T6)
- Stainless steel parts and hardware (Grade 304)
- Caster Wheels – 100mm 4 x locking
- Commode system
- Tool-free assembly and easy fold down

### SPECIFICATIONS

- Chair Width – over caster wheels 550mm
- Chair Width – over arms 550mm

- Chair Width – inside arms 450mm
- Chair Length – including footrest 840mm
- Backrest Width 448mm
- Backrest Height from Cushion to top 400mm
- Toilet Height clearance (maximum) 470mm
- Toilet Width clearance at narrowest 445mm (at casters)
- Toilet clearance can depend upon the shape of the toilet pan
- Seat Cushion Dimensions 420mm W x 420mm D
- Seat Height from Floor 530mm to 580mm (Adjustable)
- Seat Cushion Opening 175mm W x 335mm L
- Box Dimensions 780mm x 550mm x 210mm
- Box Gross Weight 9.5kg
- Actual Weight of unit 8.9kg

---

## WILL THIS FIT YOUR BATHROOM?

Use Showerbuddy's Bathcheck app to see which Showerbuddy chair will work for you - Click Here!

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)