



SHOWERBUDDY SB7E ECO TRAVELLER

\$436.00 Original price was: \$436.00.\$392.00 Current price is: \$392.00.Sale Available Online Only!



SKU: SB7e

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The ShowerBuddy SB7e Eco Traveller folding shower commode is lightweight and designed to be perfect for use in tight spaces, or for taking with you whilst travelling / storing away when not in use.

Fully height adjustable, the SB7e can roll straight over most toilets and into the shower, greatly reducing manual handling for carers, and simplifying transfers. The SB7e can also be used as a commode chair, with the bucket being removed and attached from the rear.

FEATURES

- Weight capacity 120kg
- · Rolls over most standard toilets
- Arms fixed height
- Adjustable Footrests removable and swing-away, height-adjustable
- No rust aluminium frame (Grade 6061 T6)
- Stainless steel parts and hardware (Grade 304)
- Caster Wheels 100mm 4 x locking
- Commode system
- Tool-free assembly and easy fold down

SPECIFICATIONS

- Chair Width over caster wheels 550mm
- Chair Width over arms 550mm



- Chair Width inside arms 450mm
- Chair Length including footrest 840mm
- Backrest Width 448mm
- Backrest Height from Cushion to top 400mm
- Toilet Height clearance (maximum) 470mm
- Toilet Width clearance at narrowest 445mm (at casters)
- Toilet clearance can depend upon the shape of the toilet pan
- Seat Cushion Dimensions 420mm W x 420mm D
- Seat Height from Floor 530mm to 580mm (Adjustable)
- Seat Cushion Opening 175mm W x 335mm L
- Box Dimensions 780mm x 550mm x 210mm
- Box Gross Weight 9.5kg
- Actual Weight of unit 8.9kg

WILL THIS FIT YOUR BATHROOM?

<u>Use Showerbuddy's Bathcheck app to see which Showerbuddy chair will work for you - Click Here!</u>