


SHOWERBUDDY SCP2 OPEN FRONT COMMODO CUSHION



SKU: N/A

PRODUCT DESCRIPTION

Purpose & Benefits

The **Showerbuddy SCP2 Open Front Commode Cushion** is designed to provide comfortable seating for smaller users during showering and toileting. Made from firmer PU foam with a rear incline, it supports a stable seated position while the medium-sized open front commode opening allows easier access during toileting. The cushion snaps securely into place using locating lugs and is compatible with selected Showerbuddy chairs and the P1 Petite accessory. It also comes as a standard component with the P2 Petite Buddy accessory, making it a practical solution for children and smaller users with mobility needs.

Features

- Medium-sized open front commode opening.
- Made from firmer PU foam with a rear incline.
- Snaps securely into place with locating lugs.
- Suitable for children and smaller users with mobility needs.
- Easy to clean PU foam construction.
- Compatible with selected Showerbuddy models.

- Supplied as standard with the P2 Petite Buddy accessory.

Specifications

Product Code	SCP2
Material	PU Foam with PU Cover
Commode Opening	Medium
Dimension A	14.5 cm
Dimension S	7.2 cm
Dimension T	15.5 cm
Compatible Models	SB1, SB2, SB2T, SB3T and P1 Petite accessory
Overlay	No
Insert	No

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au