

SHOWERBUDDY SEAT CUSHION WITH SMALLER COMMODE OPENING



SKU: N/A

PRODUCT DESCRIPTION

Purpose & Benefits

The **Showerbuddy Seat Cushion With Smaller Commode Opening** is designed to provide improved comfort and stability for petite users during showering and toileting. Made from firmer PU foam, the cushion features a greater incline towards the backrest to help support and stabilise the user while seated. The smaller commode opening offers a better fit for petite users, while the tool-free snap-on design allows for quick installation and removal. It is compatible with selected Showerbuddy models and is supplied as standard with the Petite accessory P1.

Features

- Designed for petite Showerbuddy users.
- Smaller commode opening for improved positioning.
- Firmer PU foam construction for added support.
- Greater rearward cushion slope enhances stability.
- Tool-free snap-on fitting with locating lugs.
- Compatible with SB1, SB2, SB2T, 3T and P2 Petite accessory.

- Supplied as standard with Petite accessory P1.

Specifications

Product Code	SCP1
Foam & Fabric	PU - PU
Dimension A	14.5 cm
Dimension S	7.2 cm
Dimension T	9 cm
Overlay	SCOC-P
Insert	No
Compatible Models	SB1, SB2, SB2T, 3T and P2 Petite accessory
Installation	Tool-free snap-on fit

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au