



SHOWERBUDDY SHOWERGLYDE MODEL 3



SKU: SG3

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **ShowerGlyde Model 3** is a compact sliding transfer bench designed to make bathing safer and easier for individuals with limited mobility. Ideal for bathrooms with vanities or tight spaces beside the tub, it features shorter rails and a swivel seat option to maximise accessibility without compromising comfort or stability. This innovative bench allows smooth, secure transfers in and out of the shower, helping maintain independence and dignity.

FEATURES

- Compact design ideal for smaller bathrooms or tubs beside vanities
- Sliding seat allows smooth transfer into and out of the tub
- Locks securely at two points – inside and outside the tub
- Rust-proof aluminium rails with stainless steel hardware
- Plastic moulded seat with multiple insert options
- Adjustable legs with non-slip feet for enhanced stability
- Swivel seat option available for improved manoeuvrability
- Includes commode seat, solid insert, and splatter guard
- Grab bar and seat belt included as standard
- Easy to reconfigure for left- or right-hand setup

SPECIFICATIONS

Specification	Detail
Overall Length	935 mm

Specification	Detail
Tub Clearance Height (Standard)	495 mm (19.5")
Tub Clearance Height (with Extenders)	Up to 635 mm (25")
Seat Options	Commode seat, solid insert, splatter guard, swivel option
Weight Capacity	150 kg
Unit Weight	10.5 kg
Colours	Black, White with Red Accents
Finish	Powder-coated aluminium rails with stainless fittings
Seat	Plastic moulded seat with optional padded cover
Armrests	Flip-back armrests (optional)
Grab Bar	Included on wall-side rail

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
 Brighton SA 5048
 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
 Ave, Elizabeth Vale SA 5112
 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
 Reynella SA 5161
 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
 03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
 03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au