



SHOWERBUDDY SHOWERGLYDE MODEL 3

\$848.00 Original price was: \$848.00.<u>\$763.20</u>Current price is: \$763.20.



SKU: SG3

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **ShowerGlyde Model 3** is a compact sliding transfer bench designed to make bathing safer and easier for individuals with limited mobility. Ideal for bathrooms with vanities or tight spaces beside the tub, it features shorter rails and a swivel seat option to maximise accessibility without compromising comfort or stability. This innovative bench allows smooth, secure transfers in and out of the shower, helping maintain independence and dignity.

FEATURES

- Compact design ideal for smaller bathrooms or tubs beside vanities
- Sliding seat allows smooth transfer into and out of the tub
- Locks securely at two points inside and outside the tub
- Rust-proof aluminium rails with stainless steel hardware
- Plastic moulded seat with multiple insert options
- Adjustable legs with non-slip feet for enhanced stability
- Swivel seat option available for improved manoeuvrability
- Includes commode seat, solid insert, and splatter guard
- Grab bar and seat belt included as standard
- Easy to reconfigure for left- or right-hand setup

SPECIFICATIONS

Specification		Detail
Overall Length	935 mm	



Specification Detail

Tub Clearance Height (Standard) 495 mm (19.5")

Tub Clearance Height (with Extenders)

Up to 635 mm (25")

Seat Options Commode seat, solid insert, splatter guard, swivel

option

Weight Capacity 150 kg
Unit Weight 10.5 kg

Colours Black, White with Red Accents

Finish Powder-coated aluminium rails with stainless fittings

Seat Plastic moulded seat with optional padded cover

Armrests Flip-back armrests (optional)

Grab Bar Included on wall-side rail