


SHOWERBUDDY SOFT FOAM CUSHION XL VINYL OVERLAY



SKU: SCSF1-3 XL

PRODUCT DESCRIPTION

Purpose & Benefits

The **Showerbuddy Soft Foam Cushion XL Vinyl Overlay** is designed to provide a softer seating surface for compatible Showerbuddy chairs while increasing the seating width for larger users. Made from medium-density foam with a machine-fitted vinyl overlay, it offers a comfortable seating experience while remaining easy to clean. The push-fit design with locating lugs allows quick installation and removal, making it a practical replacement or accessory for compatible Showerbuddy models. When paired with compatible replacement arms, it increases the seating width to accommodate a wider seating position.

Features

- Medium-density foam with a 40 mm thickness.
- Machine-fitted vinyl overlay with heat-sealed edge joins.
- Push-fit installation with locating lugs.
- Provides a softer comfort level.
- Easy-to-clean vinyl surface.
- Increases seating width to 56 cm when used as intended.

- Compatible with selected Showerbuddy models.

Specifications

Product Code	SCSF1-3 XL
Foam Type	Soft
Cover Material	Vinyl
Foam Thickness	4 cm
Seat Width	56 cm
Dimension A	18 cm
Dimension S	10 cm
Dimension T	9 cm
Overlay	No
Insert	Yes
Compatible Models	SB1, SB2, SB2T and SB3T/3T-C

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au