



## SHOWERBUDDY TRANSIT BUDDY - LIGHTWEIGHT TRAVEL CHAIR

**\$471.00**



**SKU: TB1**

## PRODUCT DESCRIPTION

### Purpose & Benefits of the ShowerBuddy Transit Buddy

The Transit Buddy TB1 is a lightweight travel chair designed to provide mobility assistance for individuals with limited mobility. Weighing just 8.4 kg, it is easily foldable, making it ideal for travel and storage. This chair facilitates smooth and comfortable mobility over short distances, enhancing accessibility and independence during outings to parks, shopping centres, or while travelling.

### Features

- **Lightweight and Foldable Design** – At just 8.4 kg, the chair folds compactly to fit into car boots, ensuring convenience for users on the go.
- **Durable Construction** – Crafted from no-rust aluminium (Grade 6061 T6), ensuring longevity and reliability.
- **Safety Mechanisms** – Equipped with a lockable handbrake system for secure stationary positioning.
- **User-Friendly Assembly** – Requires no tools for assembly and easily folds down for quick setup and storage.
- **Comfort Features** – Includes fixed-height, flip-up arms and foldable footrests to accommodate user needs.

### Specifications

- **Weight Capacity:** Supports up to 100 kg.
- **Dimensions:**
  - Chair Height: 880 mm
  - Chair Length: 925 mm
  - Width Over Arms: 540 mm

- Seat Height: 48 cm (19")
- Seat Width: 38 cm (15")
- Seat Depth: 41 cm (16")

- **Wheel Sizes:**

- Front Wheels: 6"
- Rear Wheels: 12"

- **Packaging:**

- Box Dimensions: 820 x 540 x 350 mm
- Boxed Weight: 10.6 kg

The Transit Buddy TB1 is an excellent solution for those seeking a portable and reliable mobility aid, ensuring ease of movement and comfort during daily activities.