




SHOWERBUDDY TRANSIT BUDDY - LIGHTWEIGHT TRAVEL CHAIR



SKU: TB1

PRODUCT DESCRIPTION

Purpose & Benefits of the ShowerBuddy Transit Buddy

The Transit Buddy TB1 is a lightweight travel chair designed to provide mobility assistance for individuals with limited mobility. Weighing just 8.4 kg, it is easily foldable, making it ideal for travel and storage. This chair facilitates smooth and comfortable mobility over short distances, enhancing accessibility and independence during outings to parks, shopping centres, or while travelling.

Features

- **Lightweight and Foldable Design** – At just 8.4 kg, the chair folds compactly to fit into car boots, ensuring convenience for users on the go.
- **Durable Construction** – Crafted from no-rust aluminium (Grade 6061 T6), ensuring longevity and reliability.
- **Safety Mechanisms** – Equipped with a lockable handbrake system for secure stationary positioning.
- **User-Friendly Assembly** – Requires no tools for assembly and easily folds down for quick setup and storage.
- **Comfort Features** – Includes fixed-height, flip-up arms and foldable footrests to accommodate user needs.

Specifications

- **Weight Capacity:** Supports up to 100 kg.
- **Dimensions:**
 - Chair Height: 880 mm
 - Chair Length: 925 mm
 - Width Over Arms: 540 mm
 - Seat Height: 48 cm (19")

- Seat Width: 38 cm (15")
- Seat Depth: 41 cm (16")
- **Wheel Sizes:**
 - Front Wheels: 6"
 - Rear Wheels: 12"
- **Packaging:**
 - Box Dimensions: 820 x 540 x 350 mm
 - Boxed Weight: 10.6 kg

The Transit Buddy TB1 is an excellent solution for those seeking a portable and reliable mobility aid, ensuring ease of movement and comfort during daily activities.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au