

MANGAR SIT-U-UP WITH AIRFLO 12



SKU: MPCA130500

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Sit-U-Up with Airflo 12 is a simple-to-use pillow lifter that helps to raise and lower a person in bed. It is designed to be used by a single caregiver, and it can be used to raise the user either on their back or on their side.

The Sit-U-Up is operated by low air pressure, and it is powered by the Airflo 12 compressor. The compressor is lightweight and portable, and it can be used with a variety of other Mangar products.

FEATURES

- Simple to use, even for a single caregiver
- Can be used to raise the user on their back or on their side
- Operated by low air pressure
- Powered by the Airflo 12 compressor
- Lightweight and portable
- Durable and easy to clean

The Sit-U-Up with Airflo 12 is a valuable tool for caregivers who need to help people in and out of bed. It is easy to use, safe, and effective. If you are looking for a way to make bed transfers easier, the Sit-U-Up with Airflo 12 is a great option.

SPECIFICATIONS

- Maximum user weight: 150 kilograms (24 stone)
- Length: 75 centimeters
- Width: 70 centimeters
- Thickness incl. overlay (min/max): 5 centimetres/13 centimetres
- The angle of backrest: 65° max
- Air hose length: 2 meters

- Maximum inlet pressure: 0.35 bar (5 psi)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriiffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199
03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168
03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820
1300 633 468

support@medimarnt.com.au