

MANGAR SIT-U-UP WITH AIRFLO 12

\$1,814.00



SKU: MPCA130500



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Sit-U-Up with Airflo 12 is a simple-to-use pillow lifter that helps to raise and lower a person in bed. It is designed to be used by a single caregiver, and it can be used to raise the user either on their back or on their side.

The Sit-U-Up is operated by low air pressure, and it is powered by the Airflo 12 compressor. The compressor is lightweight and portable, and it can be used with a variety of other Mangar products.

FEATURES

- Simple to use, even for a single caregiver
- Can be used to raise the user on their back or on their side
- Operated by low air pressure
- Powered by the Airflo 12 compressor
- Lightweight and portable
- Durable and easy to clean

The Sit-U-Up with Airflo 12 is a valuable tool for caregivers who need to help people in and out of bed. It is easy to use, safe, and effective. If you are looking for a way to make bed transfers easier, the Sit-U-Up with Airflo 12 is a great option.

SPECIFICATIONS

- Maximum user weight: 150 kilograms (24 stone)
- Length: 75 centimeters
- Width: 70 centimeters
- Thickness incl. overlay (min/max): 5 centimetres/13 centimetres
- The angle of backrest: 65° max

- Air hose length: 2 meters
- Maximum inlet pressure: 0.35 bar (5 psi)