



SLIDE N DINE DUAL CONTROL UTILITY CHAIR

\$1,250.00



SKU: NCM-4LW-2019C-H

PRODUCT DESCRIPTION

Purpose & Benefits

The Slide n' Dine is Australia's only dual-control dining and mobility chair. It features a unique braking system that can be released by the user with a simple lever beside their seat, or by a carer from behind the chair. This ensures maximum convenience and safety during seated transfers. Available in two frame heights, it accommodates a wide range of users. The tilt-up arms can be mounted at three height settings for personalised comfort. Choose between a standard PU foam seat or an upgraded Gelteq® seat, both with an optional coccyx cutout for pressure relief.

Features

- Dual control brake
- 4 x remote release, auto-locking castors
- Standard seat with 40mm high-resilience PU foam
- Stitch-free vinyl upholstery for easy cleaning
- Fold-up arms for easy side transfer
- Standard arm length: 340mm
- Arms adjustable to 3 height positions

- White powdercoat finish frame
- Upholstery colour options:
 - Dal Clipper Light Grey
 - Dal Clipper Light Blue
- Available in two frame heights:
 - 490mm seat height (standard)
 - 450mm seat height (low)

Options

- Extra long 390mm arms
- Gelteq™ memory foam seat upgrade
- Coccyx cutout available on both seat types
- 40mm frame height extenders
- Custom powdercoat colour options
- 21 custom upholstery options including:
 - Laine Silk Vinyl (11 colours)
 - Warwick Tritan V Healthcare Fabric (10 colours)

Specifications

- Warranty: 1 year
- Maximum user weight: 125kg

- Seat width: 450mm
- Seat depth: 440mm
- Seat height: 490mm (standard) or 450mm (low frame)