

## INVACARE SOFTFORM MAXIGLIDE MATTRESS



**SKU:** SPGM

### PRODUCT DESCRIPTION

### PURPOSE AND BENEFITS

The Softform Premier MaxiGlide Mattress features an innovative glide mechanism that significantly reduces shear and friction forces when using an electric profiling bed. The top surface provides excellent pressure redistributing properties with a high level of comfort and support, whilst the base of the mattress uniquely separates and conforms to the bedframe.

The Softform Premier MaxiGlide features a deep-contoured insert pad constructed from 'keyhole' cut castellated foam cells, which move independently to maximise the body contact area and optimise pressure reduction.

### KEY FEATURES

- Patented "Glide" Design.
- Modular Design: Easy replacement of all parts.
- Two-way stretch vapour-permeable cover: Water-proof polyurethane with welded seams.
- Fully detachable top cover: Toughened PU coated base.
- No turning required: Reduced manual handling.

### SPECIFICATIONS

<b>Brand</b>	Invacare
<b>Width</b>	Single 880mm King Single 1070mm

<b>Length</b>	1970mm
<b>Depth</b>	150mm
<b>Max SWL</b>	247.6kg
<b>Warranty</b>	4 years on cover – 8 years on foam
<b>Manufacturer's part number</b>	SPGM197/88
<b>Size</b>	Single, King Single

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
 Brighton SA 5048  
 08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
 Ave, Elizabeth Vale SA 5112  
 08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
 Reynella SA 5161  
 08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
 03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
 03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
 1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)