

## INVACARE SOFTFORM MAXIGLIDE MATTRESS

**\$1,680.00 - \$2,254.00**



**SKU:** SPMG

### PRODUCT DESCRIPTION

### PURPOSE AND BENEFITS

The Softform Premier MaxiGlide Mattress features an innovative glide mechanism that significantly reduces shear and friction forces when using an electric profiling bed. The top surface provides excellent pressure redistributing properties with a high level of comfort and support, whilst the base of the mattress uniquely separates and conforms to the bedframe.

The Softform Premier MaxiGlide features a deep-contoured insert pad constructed from 'keyhole' cut castellated foam cells, which move independently to maximise the body contact area and optimise pressure reduction.

### KEY FEATURES

- Patented "Glide" Design.
- Modular Design: Easy replacement of all parts.
- Two-way stretch vapour-permeable cover: Water-proof polyurethane with welded seams.
- Fully detachable top cover: Toughened PU coated base.
- No turning required: Reduced manual handling.

### SPECIFICATIONS

<b>Brand</b>	Invacare
<b>Width</b>	Single 880mm King Single 1070mm

<b>Length</b>	1970mm
<b>Depth</b>	150mm
<b>Max SWL</b>	247.6kg
<b>Warranty</b>	4 years on cover – 8 years on foam
<b>Manufacturer's part number</b>	SPGM197/88
<b>Size</b>	Single, King Single