



INVACARE SOFTFORM PREMIER MATTRESS

\$1,383.00 - \$1,976.00 Price range: \$1,383.00 through \$1,976.00



SKU: SPM

PRODUCT DESCRIPTION

PURPOSE AND BENEFITS

The Softform Premier ensures a sustained high degree of protection against tissue damage through the ingenious deep contoured insert pad which maximises body contact area and thereby optimises pressure reduction. The single-piece foam “U” core provides effective side-wall support to facilitate patient transfer and is suitable for use on all profiling beds. The cover fabric conforms to the deeply contoured foam castellations when the bed is profiled to minimise rucking on the top surface and reduces potential tissue damage. The toughened PU coated base prolongs the longevity of the mattress while the double-coated Poly-urethane sides reduce the risk of strike-through, hydrolysis or relamination.

KEY FUNCTIONS

- High frequency welded seams and zip - reduces risk of fluid ingress
- Full length waterfall flap conceals the zip - reduces risk of fluid ingress through zip
- No need to turn over - reduces manual handling
- Two-way stretch poly-urethane waterproof cover
- Deep contoured insert pad - enhanced pressure care area
- Single piece foam “U” core provides side wall support
- Toughened PU coated base prolongs mattress longevity

SPECIFICATIONS

- Width: 88cm - 107cm
- Height: 15cm
- Length: 197cm
- Product Weight: 14kg

- Warranty: Cover - 4 years | Foam - 8 years
- S.W.L: 247kg