




SOFTPRO FUNCTIONAL RESTING HAND SPLINT

\$270.00 - \$308.00 Price range: \$270.00 through \$308.00
Online Only



SKU: OCS-SP31FFRLL-H

PRODUCT DESCRIPTION

Purpose & Benefits of the SoftPro Functional Resting Hand Splint

The SoftPro Functional Resting Hand Splint is designed to support, protect, and correctly position the wrist, hand, and fingers for individuals experiencing mild to moderate contractures or requiring immobilisation due to conditions like rheumatoid arthritis. This splint aids in maintaining functional alignment, reducing pain, and preventing further deformity by providing gentle, consistent support.

Features

- **Bend-to-Fit Design** – Easily adjustable without the need for tools, allowing customisation to the user's specific anatomical requirements.
- **Extra-Long Wrist Strap** – Ensures the wrist is held in the proper position while applying a gentle stretch to promote optimal alignment.
- **Colour-Coded Straps** – Facilitate easy and correct application of the splint, enhancing user compliance and effectiveness.
- **Antimicrobial Material** – Constructed with materials that help prevent odour and infection, ensuring hygiene and comfort during extended use.
- **Thumb Support** – Features a bendable component to support, position, and correct the alignment of the thumb as needed.

Specifications

Size	Width Of MP Joints	Wrist To Fingertip Length
Small	Up to 7.6cm	Up to 15.2cm
Medium	7.6-8.9cm	15.2-20.3cm

Large

8.9-10.2cm

20.3-24.1cm

To determine the correct size, measure the width across the metacarpophalangeal (MP) joints and the length from the wrist to the tip of the middle finger.

The SoftPro Functional Resting Hand Splint offers an effective solution for individuals seeking to maintain hand functionality, alleviate pain, and prevent further contractures through gentle support and proper positioning.