



### SOFTPRO FUNCTIONAL RESTING HAND SPLINT



SKU: OCS-SP31FFRLL-H

### PRODUCT DESCRIPTION

# Purpose & Benefits of the SoftPro Functional Resting Hand Splint

The SoftPro Functional Resting Hand Splint is designed to support, protect, and correctly position the wrist, hand, and fingers for individuals experiencing mild to moderate contractures or requiring immobilisation due to conditions like rheumatoid arthritis. This splint aids in maintaining functional alignment, reducing pain, and preventing further deformity by providing gentle, consistent support.

### **Features**

- **Bend-to-Fit Design** Easily adjustable without the need for tools, allowing customisation to the user's specific anatomical requirements.
- Extra-Long Wrist Strap Ensures the wrist is held in the proper position while applying a gentle stretch to promote optimal alignment.
- **Colour-Coded Straps** Facilitate easy and correct application of the splint, enhancing user compliance and effectiveness.
- **Antimicrobial Material** Constructed with materials that help prevent odour and infection, ensuring hygiene and comfort during extended use.
- **Thumb Support** Features a bendable component to support, position, and correct the alignment of the thumb as needed.

# **Specifications**

SizeWidth Of MP JointsWrist To Fingertip LengthSmallUp to 7.6cmUp to 15.2cmMedium7.6-8.9cm15.2-20.3cm



Large

8.9-10.2cm

20.3-24.1cm

To determine the correct size, measure the width across the metacarpophalangeal (MP) joints and the length from the wrist to the tip of the middle finger.

The SoftPro Functional Resting Hand Splint offers an effective solution for individuals seeking to maintain hand functionality, alleviate pain, and prevent further contractures through gentle support and proper positioning.

For further information go to medimart.com.au or give us a call: 1300 633 063

### SA

# **Brighton Showroom**

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

# **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

# **Reynella Showroom**

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

# **VIC**

#### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451

frankston@medimart.com.au

#### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au

# NT

#### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 <a href="mailto:support@medimarnt.com.au">support@medimarnt.com.au</a>