

## SOVEREIGN S3 PRESSURE CARE MATTRESS

**\$1,447.00 - \$1,452.00** Price range: \$1,447.00 through \$1,452.00



**SKU:** S3SQL1000-H



## PRODUCT DESCRIPTION

### Purpose & Benefits of the Sovereign S3 Pressure Care Mattress

The Sovereign S3 Pressure Care Mattress is a premium solution for individuals at high risk of pressure injuries, offering exceptional support, durability, and comfort. Engineered with advanced foam technology, the mattress provides optimal pressure redistribution while remaining lightweight for easy handling. Ideal for homecare, aged care, and clinical settings, the Sovereign S3 ensures long-term comfort and effective pressure care for users up to 250 kg.

### Features

- **Triple-Layer Foam Construction:** Combines three high-quality foam layers to deliver superior pressure relief, postural support, and stability for high-risk users.
- **Lightweight Design:** Despite its robust 250 kg weight capacity, the mattress remains lightweight, making it easy for caregivers to handle during repositioning and transport.
- **Firm Edge Support:** Reinforced side walls enhance stability during transfers, reduce the risk of roll-out, and maintain mattress structure.
- **Temperature Regulation:** Incorporates breathable materials to promote airflow, ensuring a cool and comfortable sleep surface.
- **Vapour-Permeable, Waterproof Cover:** The multi-stretch cover is waterproof, breathable, and easy to clean, supporting infection control and maintaining hygiene.
- **Durable and Long-Lasting:** Built to withstand consistent use while maintaining its integrity, ensuring reliability in high-demand environments.

### Specifications

- **Weight Capacity:** 250 kg

- **Mattress Dimensions:** Available in Single, King Single, Double, and custom sizes
- **Foam Layers:** Triple-layer pressure care foam
- **Risk Level:** Suitable for users at high risk of pressure injuries
- **Cover Material:** Vapour-permeable, waterproof, and multi-stretch fabric
- **Cleaning:** Wipe clean with a damp cloth or disinfectant
- **Usage:** Suitable for homecare, aged care, and clinical environments

The Sovereign S3 Pressure Care Mattress delivers advanced pressure care, stability, and comfort, providing a reliable solution for individuals requiring long-term support in high-risk care scenarios.