

SOVEREIGN S7 PRESSURE CARE MATTRESS





SKU: S7D1000PU-H

PRODUCT DESCRIPTION

Purpose & Benefits of the Sovereign S7 Pressure Care Mattress

The Sovereign S7 Pressure Care Mattress is a premium solution designed for individuals at very high risk of pressure injuries. Featuring advanced triple-layer foam technology and reinforced support, the S7 delivers exceptional comfort, durability, and stability for users requiring intensive care. With a high weight capacity of up to 250 kg, this mattress is ideal for homecare, aged care, and clinical environments, ensuring user safety and long-term performance.

Features

- **Triple-Layer Foam System**: Combines high-density support foam, a comfort mid-layer, and a visco-elastic memory foam top layer to provide superior pressure redistribution and postural support.
- **High Weight Capacity**: Engineered to support users up to 250 kg, making it ideal for bariatric and high-dependency care needs.
- **Reinforced Edge Support**: The firm side walls ensure stability during transfers, reduce the risk of roll-out, and maintain the mattress structure over time.
- **Visco-Elastic Memory Foam**: The top layer conforms to the user's body, distributing weight evenly and minimising pressure points for optimal comfort.
- **Temperature Regulation**: Open-cell memory foam promotes airflow to regulate temperature and prevent heat build-up, keeping the user cool and comfortable.
- **Vapour-Permeable, Waterproof Cover**: The multi-stretch, waterproof cover is easy to clean and enhances infection control, ensuring a hygienic sleep surface.
- **Built for Longevity**: Designed for long-term use, the S7 offers reliable performance and durability under high-demand conditions.



Specifications

- Weight Capacity: 250 kg
- Mattress Dimensions: Available in Long Single, Long King Single, Double, and Long Double
- Foam Layers: Triple-layer with visco-elastic memory foam
- Risk Level: Suitable for individuals at very high risk of pressure injuries
- Cover Material: Waterproof, vapour-permeable, multi-stretch fabric
- Cleaning: Wipe clean with a damp cloth or disinfectant
- **Usage**: Ideal for bariatric care, aged care, and clinical environments requiring advanced pressure care

The Sovereign S7 Pressure Care Mattress offers a robust and effective solution for individuals requiring intensive care, combining comfort, stability, and durability to enhance user safety and well-being.

For further information go to medimart.com.au or give us a call: 1300 633 063

SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 319903 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au

Page: 2

