

SOVEREIGN S7 PRESSURE CARE MATTRESS

\$2,243.00 - \$4,125.00 Price range: \$2,243.00 through \$4,125.00



SKU: S7D1000PU-H



PRODUCT DESCRIPTION

Purpose & Benefits of the Sovereign S7 Pressure Care Mattress

The Sovereign S7 Pressure Care Mattress is a premium solution designed for individuals at very high risk of pressure injuries. Featuring advanced triple-layer foam technology and reinforced support, the S7 delivers exceptional comfort, durability, and stability for users requiring intensive care. With a high weight capacity of up to 250 kg, this mattress is ideal for homecare, aged care, and clinical environments, ensuring user safety and long-term performance.

Features

- **Triple-Layer Foam System:** Combines high-density support foam, a comfort mid-layer, and a visco-elastic memory foam top layer to provide superior pressure redistribution and postural support.
- **High Weight Capacity:** Engineered to support users up to 250 kg, making it ideal for bariatric and high-dependency care needs.
- **Reinforced Edge Support:** The firm side walls ensure stability during transfers, reduce the risk of roll-out, and maintain the mattress structure over time.
- **Visco-Elastic Memory Foam:** The top layer conforms to the user's body, distributing weight evenly and minimising pressure points for optimal comfort.
- **Temperature Regulation:** Open-cell memory foam promotes airflow to regulate temperature and prevent heat build-up, keeping the user cool and comfortable.
- **Vapour-Permeable, Waterproof Cover:** The multi-stretch, waterproof cover is easy to clean and enhances infection control, ensuring a hygienic sleep surface.
- **Built for Longevity:** Designed for long-term use, the S7 offers reliable performance and durability under high-demand conditions.

Specifications

- **Weight Capacity:** 250 kg
- **Mattress Dimensions:** Available in Long Single, Long King Single, Double, and Long Double
- **Foam Layers:** Triple-layer with visco-elastic memory foam
- **Risk Level:** Suitable for individuals at very high risk of pressure injuries
- **Cover Material:** Waterproof, vapour-permeable, multi-stretch fabric
- **Cleaning:** Wipe clean with a damp cloth or disinfectant
- **Usage:** Ideal for bariatric care, aged care, and clinical environments requiring advanced pressure care

The Sovereign S7 Pressure Care Mattress offers a robust and effective solution for individuals requiring intensive care, combining comfort, stability, and durability to enhance user safety and well-being.