

## STAEBEL BARIATRIC SHOWER TROLLEY





**SKU:** 12165

### PRODUCT DESCRIPTION

# **Purpose & Benefits**

The **Staebel Bariatric Shower Trolley** is a high-quality, heavy-duty solution designed for safe and comfortable hygiene care of bariatric patients. With a Safe Working Load of 250 kg, this shower trolley provides stability, safety and ease of use in hospitals, aged care facilities and clinical environments. It combines strength, comfort and functionality to ensure optimal patient support and efficient caregiver operation during showering and personal care routines.

#### **Features**

- Spacious stretcher top with 90 cm internal width and 107 cm external width for patient comfort and accessibility
- Hydraulic height adjustment from 54 cm to 97 cm for smooth transfers and ergonomic caregiver
  use
- Large 150 mm castors for easy manoeuvrability across various floor surfaces
- Central braking system with directional wheel lock for improved safety and control during transfers
- Foldable powder-coated stainless steel side rails with safety release mechanism for secure patient handling



- Water-resistant 27 cm thick shower mattress with integrated drainage and flexible hose (645 mm to 2000 mm)
- · Semi-cylinder head pillow for additional head and neck support during showering
- Durable metal base cover for protection and longevity in wet environments
- Stable and secure design, ideal for routine hygiene and rehabilitation care

# **Specifications**

Specification	Details
---------------	---------

Model Staebel Bariatric Shower Trolley

SKU 1216ST Safe Working Load (SWL) 250 kg

Height Adjustment Range 54 cm - 97 cm (hydraulic operation)

Internal Width 90 cm External Width 107 cm

Mattress Thickness 27 cm water-resistant mattress with integrated drainage

Hose Length 645 mm - 2000 mm flexible drainage hose

Side Rails 18.5 cm high, foldable, powder-coated stainless steel with safety

release

Castors 150 mm heavy-duty castors with central braking and steering lock

Base Cover Durable metal cover for protection and longevity

Trolley Weight 110 kg