

STAEBEL BARIATRIC SHOWER TROLLEY



SKU: 1216ST

PRODUCT DESCRIPTION

Purpose & Benefits

The **Staebel Bariatric Shower Trolley** is a high-quality, heavy-duty solution designed for safe and comfortable hygiene care of bariatric patients. With a Safe Working Load of 250 kg, this shower trolley provides stability, safety and ease of use in hospitals, aged care facilities and clinical environments. It combines strength, comfort and functionality to ensure optimal patient support and efficient caregiver operation during showering and personal care routines.

Features

- Spacious stretcher top with 90 cm internal width and 107 cm external width for patient comfort and accessibility
- Hydraulic height adjustment from 54 cm to 97 cm for smooth transfers and ergonomic caregiver use
- Large 150 mm castors for easy manoeuvrability across various floor surfaces
- Central braking system with directional wheel lock for improved safety and control during transfers
- Foldable powder-coated stainless steel side rails with safety release mechanism for secure patient handling

- Water-resistant 27 cm thick shower mattress with integrated drainage and flexible hose (645 mm to 2000 mm)
- Semi-cylinder head pillow for additional head and neck support during showering
- Durable metal base cover for protection and longevity in wet environments
- Stable and secure design, ideal for routine hygiene and rehabilitation care

Specifications

Specification	Details
Model	Staebel Bariatric Shower Trolley
SKU	1216ST
Safe Working Load (SWL)	250 kg
Height Adjustment Range	54 cm – 97 cm (hydraulic operation)
Internal Width	90 cm
External Width	107 cm
Mattress Thickness	27 cm water-resistant mattress with integrated drainage
Hose Length	645 mm – 2000 mm flexible drainage hose
Side Rails	18.5 cm high, foldable, powder-coated stainless steel with safety release
Castors	150 mm heavy-duty castors with central braking and steering lock
Base Cover	Durable metal cover for protection and longevity
Trolley Weight	110 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North

Brighton SA 5048

08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice

Ave, Elizabeth Vale SA 5112

08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,

Reynella SA 5161

08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au