


## STAEBEL JUNIOR SHOWER TROLLEY



**SKU:** 1214ST



## PRODUCT DESCRIPTION

## PURPOSE & BENEFITS OF THE STAEBEL JUNIOR SHOWER TROLLEY

The Staebel Junior Shower Trolley is a robust and reliable solution designed to assist caregivers in providing safe and comfortable showering for individuals with limited mobility. Ideal for use in hospitals, care facilities, and home care environments, this trolley allows for easy patient transfers and bathing while maintaining dignity and safety. The height-adjustable design and durable construction provide comfort for both the patient and caregiver, ensuring a smooth and efficient showering process.

## FEATURES

- **Height Adjustable:** The trolley can be adjusted in height, allowing caregivers to set the optimal working level, reducing strain during patient transfers and showering.
- **Durable Construction:** Built with a sturdy frame and high-quality materials, the Staebel Junior can support users up to 170 kg, making it suitable for a wide range of individuals.
- **Easy Patient Transfer:** The trolley's design allows for smooth and safe transfers, minimising the risk of injury to both the patient and caregiver during the bathing process.
- **Waterproof and Easy to Clean:** Constructed from waterproof materials, the trolley is easy to clean and maintain, ensuring hygiene in clinical and home care settings.
- **Comfortable and Safe Design:** Equipped with padded surfaces for patient comfort and side rails for added safety during transfers and bathing.

## SPECIFICATIONS

- Shower trolley with one-section trolley top
- Trolley top length: 150 cm, trolley top width: 60 cm

- Maximum dimensions: 167 cm (length) and 76 cm (width)
- Adjustable 5 trolley top inclination
- Hydraulic lowering/raising of trolley top via a console from 54cm to 97cm
- Castors with diameter of 150 mm
- Central braking system
- Easy steering system (for ease of maneuvering trolleys front left castor can be locked in direction of steering)
- Side rails foldable, height: 18,5 cm with safety release mechanism (material - powder coated stainless steel)
- Metal base cover
- Water-resistant shower mattress (height: 27 cm)
- Mattress integrated drain from 645-2000 mm long flexible hose
- Water-resistant head pillow semi-cylinder
- Trolley weight: 79 kg
- Safe working load: 170 kg

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

**Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)