



STAEBEL TOTAL M9500 BARIATRIC FOAM MATTRESS - DOUBLE

\$2,574.00



SKU: M9500D

PRODUCT DESCRIPTION

Purpose & Benefits

The **Staebel Total M9500 Bariatric Foam Mattress - Double** provides advanced pressure care, comfort, and support for bariatric patients weighing up to 500 kg. Designed with a 4-layered graded pressure care core and premium foam materials, it delivers long-lasting durability and clinical-grade performance for superior patient outcomes.

Features

- 4-layered graded pressure care core for optimal comfort and support
- High-specification foam construction ensures performance and longevity
- Premium visco-grade upper immersion layer provides low-shear surface and effective pressure relief
- Secondary postural support layer maintains patient alignment and stability in various positions
- Layered firm base prevents bottoming out and enhances support during repositioning
- Triple-section, split base design optimises profiling with adjustable bed frames
- Anti-microbial treatment improves hygiene and mattress longevity

- Foam conforms to AS/NZS 4088.1 (Int):1993, FMVSS302, and SAE T2002 standards

Specifications

- **Model:** Total M9500 Bariatric Mattress
- **Size:** Double – 1960 mm (L) x 1350 mm (W) x 175 mm (H)
- **Safe Working Load (SWL):** 500 kg
- **Core Structure:** 4-layer graded pressure care foam
- **Base Type:** Triple-section split base
- **Material Treatment:** Anti-microbial foam for improved hygiene