


STANDARD FACE MASKS - 50 PACK

\$25.00



SKU: 195663



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Our masks are made of environmentally friendly and breathable material to keep you feeling fresh while safely protected from germs, dust, pollen and more. The outer non-woven fabric hydrophobic layer repels water, blood and other bodily fluids. The middle melt-down fabric filter layer stops harmful particles from entering or exiting the mask. And the inner non-woven fabric hydrophilic layer is soft and absorbent, to keep the mask dry, skin-friendly, allergy-free and comfortable to wear. The concealed metal nose strip moulds seamlessly to conform to your facial contours, so that the mask can fit various face shapes and provide better protection. The elastic ear loops are strong and durable but not tight, allowing for long-term wear without hurting your ears. Our high-quality disposable mask is not only suitable for use by the general public but also in specialised areas such as hospitals, nail salons, factories and garages. Breathe confidently when you are out and about, knowing you can rely on our disposable mask for daily protection.

FEATURES

- Disposable facial mask
- 100% brand new and high quality
- Environmentally friendly fabric
- 3-layer protection
- Anti PM2.5, fine particulate matter
- Dust and fluid-resistant
- Filters out at least 95% of particles
- Skin-friendly
- Breathable

- Easy to wear
- Adjustable metal nose strip
- Elastic ear loops
- Perfect edging
- Perfect fitting
- CE-certified