




## STANDER ASSIST-A-TRAY



SKU: 41012

## PRODUCT DESCRIPTION

### PURPOSE & BENEFITS

The Stander Assist-A-Tray is a versatile and innovative support device designed for individuals facing mobility or balance challenges. It serves as a combination of a standing aid and a convenient tray table, ideal for use with a chair, sofa, or bed.

### BENEFITS

- **Dual Functionality:** Incorporates a sturdy, adjustable tray table and a handle for support, offering a convenient surface for various activities and aid in standing.
- **Ergonomic Safety Handle:** Provides a safe and secure grip, aiding in safely sitting down or standing up, thus reducing the risk of falls.
- **Adjustable Height and Angle:** The tray table can be tailored in height and angle, accommodating activities like eating, reading, or using a laptop.
- **Swivel Tray:** The tray is capable of a 360-degree swivel for easy access and positioning.
- **Sturdy and Stable Design:** Constructed with a durable frame to support a significant amount of weight, ensuring stability.
- **Easy Assembly and Compatibility:** Designed for straightforward installation and compatible with most chairs, sofas, and beds.

The Stander Assist-A-Tray is a key aid for enhancing safety, independence, and comfort in the living space. It is especially beneficial for those with limited mobility or balance issues, making daily tasks more accessible and manageable.

## SPECIFICATIONS

- Height Adjustment of Tray: 66.04 cm - 81.28 cm
- Height Adjustment of Handle: 86.36 cm - 101.6 cm
- Handle Size: 15.24 cm H x 15.24 cm W x 3.175 cm D
- Tray Size: 40.64 cm x 45.72 cm
- Base Length Adjustment: 50.8 cm - 91.44 cm
- Base Pads: 10.16 cm x 17.78 cm
- Weight: 11.34 kg
- Weight Capacity of Handle: 113.4 kg
- Weight Capacity of Tray: 13.61 kg

## INSTRUCTION MANUAL

[Stander Assist-A-Tray Instruction Manual](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### Brighton Showroom

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### Reynella Showroom

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### Monash Showroom

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)

