



STANDER SECURITY POLE- CURVED GRAB BAR (WHITE)

\$492.00



SKU: 47003

PRODUCT DESCRIPTION

Purpose & Benefits of Stander Security Pole - Curved Grab Bar (White)

The Stander Security Pole with a Curved Grab Bar in white offers a versatile and secure solution to aid mobility and stability around the home, particularly in the bathroom. This innovative product is designed to assist individuals with limited mobility or those who require extra support when standing up from a sitting position. It's an excellent aid for seniors or anyone needing additional stability, providing a safe means to navigate daily activities independently and confidently.

Features

- **Versatile Functionality:** The pole extends from the floor to the ceiling, offering sturdy support wherever it's installed—be it next to a bed, bathtub, toilet, or chair.
- **Curved Grab Bar:** The grab bar at the top provides a multi-position handgrip that allows users to gain leverage when standing up, adding comfort and safety. It can be easily rotated and locked at multiple angles as per the user's needs.
- **Adjustable Height:** The security pole can be adjusted to fit ceiling heights from 213 to 305 cm, making it a flexible option for different room sizes.
- **No Wall Mounts Needed:** Installation does not require wall mounting, making it an excellent choice for areas where permanent fixtures are not possible. The tension-mounted design ensures easy installation and removal without damage to surfaces.
- **High Weight Capacity:** Supports up to 136 kg, providing a reliable and durable support structure.

Specifications

Specification	Value
Fits Ceiling Heights	213 - 305 cm
Weight Capacity	136 kg
Weight	9 kg
Height of Curved Grab Bar	81 - 132 cm from floor
Length of Curved Grab Bar	49.5 cm
Width of Curved Grab Bar	21.6 cm
Length of Top I-Beam	56 cm
Width of Top I-Beam	23 cm
Diameter of Circular Base	18.4 cm
Diameter of Pole	3.8 cm

Video

The Stander Security Pole – Curved Grab Bar is an indispensable tool for enhancing safety and mobility within the home, offering peace of mind for users and their families. Its robust construction and thoughtful design ensure it is a practical addition to any living space requiring added support.