



STANDEZY BED POLE



SKU: N/A

PRODUCT DESCRIPTION

Purpose & Benefits of the Standezy Bed Pole

The Standezy Bed Pole is a versatile and ergonomic mobility aid designed to provide essential support for individuals who need assistance with repositioning in bed. Perfect for use in home settings or care facilities, this bed pole offers a reliable solution for enhancing independence and safety during daily routines.

Key Features

- **Ergonomic Design:** The Standezy Bed Pole features an ergonomic handle that provides a comfortable and secure grip, allowing users to pull themselves up, turn over, or reposition with ease.
- **Adjustable Height:** The pole's height can be easily adjusted to suit the user's specific needs, ensuring optimal support and comfort during use.
- **Sturdy Construction:** Built with high-quality, durable materials, the Standezy Bed Pole is designed to withstand regular use while providing stable and reliable support.
- **Easy Installation:** The bed pole can be quickly and securely installed between the mattress and bed base, making it a convenient and practical solution for most bed types.
- **Compact and Portable:** Lightweight and easy to handle, the Standezy Bed Pole is also portable, making it an excellent choice for those who may need support in different locations.

Specifications

- **Overall Dimensions Depth:** 45 mm
- **Height:** 690 mm
- **Handle Height:** 540 mm

- **Product Weight:** 2.1 kg
- **Safe Working Load:** 300 kg

The Standezy Bed Pole is an ideal mobility aid for those seeking additional support to enhance their independence in bed. Its ergonomic design, adjustable height, and sturdy construction make it a reliable and practical choice for improving daily comfort and safety.