


## STEADY STEP® SOL STEP



**SKU:** SP-50040-H

## PRODUCT DESCRIPTION

### Purpose & Benefits of Steady Step® Sol Step

The Steady Step® Sol Step is designed to provide relief and support for individuals suffering from heel pain, plantar fasciitis, or heel spurs. This specialized orthotic device features a unique design that cushions and supports the heel, reducing impact and stress on the painful area while walking or standing. Its shock-absorbing properties help mitigate pain and discomfort, allowing for increased mobility and a quicker return to daily activities.

### Features

- **Gel Cushioning:** Incorporates a gel pad in the heel area that absorbs impact and distributes pressure evenly to alleviate pain.
- **Compression Therapy:** Offers gentle compression that improves blood flow to the heel, helping reduce swelling and providing pain relief.
- **Adjustable Straps:** Comes with adjustable straps that ensure a snug fit and prevent the device from shifting during movement.
- **Neoprene Construction:** Made from durable neoprene that is soft against the skin and provides sufficient stretch to accommodate different foot sizes comfortably.
- **Lightweight Design:** The device is lightweight and thin enough to wear under socks and in most shoes without altering the fit.
- **Ease of Use:** Simple to put on and take off, the Sol Step can be used throughout the day during various activities for continuous support.

### Specifications

- **Material:** High-quality neoprene with integrated gel pads.

- **Size:** Small, Medium, Large, Extra Large.
- **Care Instructions:** Easy to maintain; can be hand washed with mild detergent and air dried.

The Steady Step® Sol Step is ideal for those looking to alleviate heel pain associated with conditions like plantar fasciitis and heel spurs. Its therapeutic design not only offers immediate pain relief but also supports long-term healing by reducing stress on the heel during daily activities.

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

## SA

### **Brighton Showroom**

343 Brighton Rd, North  
Brighton SA 5048  
08 7079 8773

[brighton@medimart.com.au](mailto:brighton@medimart.com.au)

### **Elizabeth Vale Showroom**

Playford Health Hub, 44 John Rice  
Ave, Elizabeth Vale SA 5112  
08 7079 8777

[elizabethvale@medimart.com.au](mailto:elizabethvale@medimart.com.au)

### **Reynella Showroom**

Shop 1/141 Sherriffs Rd,  
Reynella SA 5161  
08 7079 8771

[reynella@medimart.com.au](mailto:reynella@medimart.com.au)

## VIC

### **Frankston Showroom**

111 Cranbourne Rd, Frankston VIC 3199  
03 8763 3451

[frankston@medimart.com.au](mailto:frankston@medimart.com.au)

### **Monash Showroom**

23 Winterton Road, Clayton VIC 3168  
03 7503 8354

[monash@medimart.com.au](mailto:monash@medimart.com.au)

## NT

### **Darwin Showroom**

Unit 4/524 Stuart Hwy, Winnellie NT 0820  
1300 633 468

[support@medimarnt.com.au](mailto:support@medimarnt.com.au)