



SUPA SCOOTA SUPALITE

\$3,900.00 Original price was:

~~\$3,900.00~~ \$3,510.00 Current price is: **\$3,510.00.**



SKU: SupaLite3

PRODUCT DESCRIPTION

Purpose & Benefits of the Supa Scoota Supalite

The Supa Scoota Supalite is a highly portable mobility scooter designed to offer the utmost convenience and ease of use for individuals with mobility limitations. As one of the lightest scooters on the market, the Supalite is perfect for those who travel frequently or need a compact solution for daily mobility. Its lightweight design does not compromise performance or durability, making it an excellent choice for navigating both indoor and outdoor environments.

Features

- **Ultra-Lightweight Design:** Weighing just around 20.4 kg (45 lbs) including the battery, the Supalite is exceptionally easy to lift and transport, ideal for users who frequently transition between travel and everyday use.
- **Compact and Foldable:** Features a simple folding mechanism that allows for quick and easy storage in car trunks, closets, or while travelling.
- **Adjustable Components:** Includes an adjustable tiller (steering column) and a comfortable seat that caters to the comfort and needs of the user.
- **Durable Construction:** Built with high-quality materials that ensure reliability and durability without adding unnecessary weight.
- **Efficient Battery System:** Equipped with a lithium battery that provides ample range for daily activities, ensuring users can go about their day without concerns about battery life.
- **Easy Maneuverability:** Designed for excellent maneuverability in tight spaces, making it suitable for both indoor and outdoor environments.

Specifications

- **Material:** Constructed from high-grade aluminium and other lightweight materials for durability and ease of transport.
- **Weight:** Approximately 20.4 kg (45 lbs), including the battery.
- **Battery Type:** Lithium-ion battery that offers a balance between weight and performance.
- **Dimensions:**
 - Length 970 mm
 - Width 570 mm
 - Pneumatic Tyres 200 x 65 mm
- **Safe Working Load:** Supports users up to about 120 kg.
- **Travel Range:** 10 – 13 km (Estimate only. Will vary according to terrain, user weight, battery condition, etc.)
- **Care Instructions:** Regular maintenance includes battery care, cleaning, and inspection of moving parts for wear and tear.

The Supa Scoota Supalite is ideal for users seeking a lightweight, portable mobility scooter that does not sacrifice performance for portability. Its ease of use, combined with its compact and lightweight design, makes it a preferred choice for those needing mobility assistance while maintaining an active and independent lifestyle.