




TA INDOOR WAVE POWER CHAIR



SKU: TAIWPC

PRODUCT DESCRIPTION

The TA Indoor Power Chair was developed and manufactured in Denmark by TA Service. It is a Class A Indoor/Outdoor Chair tested and approved at the Swedish Institute of Assistive Technology. It is primarily used for indoor and outdoor use around the home. The most common specification includes electric elevate and tilt, for optimum functionality.

The TA Indoor is designed to be as small and compact as possible. Being only 56cm wide (excluding seat) and 69.5cm long (excluding leg rests) it will easily negotiate doorways, narrow passageways and tighter corners with a small turning radius of only 45cm. Due to its small compact design and exceptional manoeuvrability in some situations the TA Indoor will save the inconvenience and cost of modifying doorways, bathrooms and kitchens. Being a modular system, with highly adjustable components, ensures fewer specific parts are required to reconfigure for other users specific needs.

Specifications

- Maximum User Weight: 125kg
- Lifting Height: 32cm
- Tilt: 34 degrees
- Width excl. seat: 56cm
- Length excl. seat: 69.5cm