



TA INDOOR WAVE POWER CHAIR



SKU: TAIWPO

PRODUCT DESCRIPTION

The TA Indoor Power Chair was developed and manufactured in Denmark by TA Service. It is a Class A Indoor/Outdoor Chair tested and approved at the Swedish Institute of Assistive Technology. It is primarily used for indoor and outdoor use around the home. The most common specification includes electric elevate and tilt, for optimum functionality.

The TA Indoor is designed to be as small and compact as possible. Being only 56cm wide (excluding seat) and 69.5cm long (excluding leg rests) it will easily negotiate doorways, narrow passageways and tighter corners with a small turning radius of only 45cm. Due to its small compact design and exceptional manouevrability in some situations the TA Indoor will save the inconvenience and cost of modifying doorways, bathrooms and kitchens.Being a modular system, with highly adjustable components, ensures fewer specific parts are required to reconfigure for other users specific needs.

Specifications

Maximum User Weight: 125kg

Lifting Height: 32cmTilt: 34 degrees

Width excl. seat: 56cmLength excl. seat: 69.5cm

For further information go to medimart.com.au or give us a call: 1300 633 063



SA

Brighton Showroom

343 Brighton Rd, North Brighton SA 5048 08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice Ave, Elizabeth Vale SA 5112 08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd, Reynella SA 5161 08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199 03 8763 3451 frankston@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820 1300 633 468 support@medimarnt.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168 03 7503 8354 monash@medimart.com.au