

TAB TIMER RAINBOW CONNECTION PILL BOX

~~\$77.00~~ Original price was: \$77.00. \$69.30 Current price is: \$69.30.



SKU: TT4-15

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The **Tab Timer Rainbow Connection Pill Box** is a practical and flexible medication organiser designed to simplify daily pill management. With up to four alarms per day and up to fifteen interlocking compartments, it helps users organise multiple medications, vitamins, or supplements with precision. Ideal for home, travel, or aged care, it features clear sound and visual reminders to ensure doses are never missed, supporting routine, independence, and peace of mind.

FEATURES

- Connect up to 15 compartments for customised storage and easy daily organisation
- Set up to 4 alarms per day with clear 1-minute sound alerts
- Repeating reminder for missed doses at 5-minute intervals (3 times)
- LCD display with 12/24-hour format for clear time visibility
- Made from food-grade, BPA-free Tritan plastic for safe use
- Includes 2 × AAA batteries for immediate setup
- Compact, lightweight design ideal for home use or travel
- Backed by a 12-month limited warranty (excluding batteries)

SPECIFICATIONS

- **Brand:** Tab Timer
- **Product Name:** Rainbow Connection Pill Box
- **Alarms:** Up to 4 per day
- **Compartments:** 1-15 connectable sections
- **Alarm Duration:** 1 minute, repeats 3 times after 5 minutes if missed

- **Sound Level:** 45-50 dB
- **Display:** LCD (12/24-hour format)
- **Material:** Food-grade Tritan, BPA-free
- **Weight:** 83 g
- **Power:** 2 × AAA batteries (included)
- **Timer Dimensions:** 70 mm × 70 mm × 25 mm
- **Timer with 7 Compartments Connected:** 70 mm × 300 mm × 25 mm
- **Compartment Size:** 25 mm × 25 mm × 22 mm (fits 6-8 large tablets or 10-25 small pills)
- **Warranty:** 12-month limited (excludes battery)