


TACTEE CUTLERY SYSTEM KIT



SKU: TAC-AT-KIT



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Tactee system is a set of assistive devices designed to restore autonomy and independence in the daily activities of people with mono or bilateral functional grasping deficits that aim to improve their quality of life. The Tactee system returns functional autonomy and independence in different daily activities such as eating, drinking, writing, etc.; eliminating the caregiver who's been considered indispensable so far.

When the grasping element is approached by a Tactee, tool, the magnetic force attracts it and allows a firm grip without using fingers. At the same time, its special design allows an easy release, permitting interchanging tools without any effort nor assistance, in total autonomy. The Tactee method involves the use of the system along the rehabilitation path of the users and promotes it in every context of everyday life, in an easy, decent and independent way. Tactee is a patented system based on an innovative magnetic coupling solution, between a grasping element worn by the hand's user and various interchangeable tools that guarantee high performance in autonomy and safety.

FEATURES

- **Ergonomic Design:** Handles are crafted to fit comfortably in the hand, reducing strain and enhancing grip, even for those with limited dexterity.
- **Modular Attachments:** A selection of interchangeable utensil heads, including spoons, forks, and knives, allows users to customize the cutlery to their specific needs.
- **Secure Connection:** The system features a simple yet secure mechanism for attaching and detaching the heads, ensuring they stay in place during use.

- **Durable Materials:** Constructed with high-quality, dishwasher-safe materials, the cutlery is easy to clean and designed to withstand daily use.
- **Non-Slip Handles:** The handles are textured or coated with a non-slip material to provide a stable grip, even for users with compromised hand strength.

TACTEE KIT CONTENTS

- **Grasping Element** - the flexible strap made from polyurethane that makes it both hygienic and easy to close the grasping element using the mouth. To secure the grasping element further there is a thumb strap that prevents the adapter from sliding off the hand once loaded. Can be used on the right or left hand. Available in S, M, or L.
- **Fork, Spoon & Knife** - Injection moulded and made of stainless steel and polycarbonate. The angle of the fork, spoon & knife can be adjusted by using a bending iron.
- **Universal Adapter** - can be used with a toothbrush, razor, pen or other items.
- **Bottle Holder** - has a velcro strap attached which makes it easy to attach a bottle or cup.
- **Pads** - a set of dorsal pads that can be applied on the Grasping Element for users who are pressure-sensitive.
- **Black Storage Case** - allows you to store, protect and transport the Tactee Kit contents safely.

SIZING

Available in 3 sizes. Measurement to be taken around MCP joints.

BROCHURE

[View the Medimart Tactee brochure.](#)

For further information go to **medimart.com.au** or give us a call: **1300 633 063**

SA

Brighton Showroom

343 Brighton Rd, North
Brighton SA 5048
08 7079 8773

brighton@medimart.com.au

Elizabeth Vale Showroom

Playford Health Hub, 44 John Rice
Ave, Elizabeth Vale SA 5112
08 7079 8777

elizabethvale@medimart.com.au

Reynella Showroom

Shop 1/141 Sherriffs Rd,
Reynella SA 5161
08 7079 8771

reynella@medimart.com.au

VIC

Frankston Showroom

111 Cranbourne Rd, Frankston VIC 3199

03 8763 3451

frankston@medimart.com.au

Monash Showroom

23 Winterton Road, Clayton VIC 3168

03 7503 8354

monash@medimart.com.au

NT

Darwin Showroom

Unit 4/524 Stuart Hwy, Winnellie NT 0820

1300 633 468

support@medimarnt.com.au