

TACTEE CUTLERY SYSTEM KIT



\$592.00 Original price was: \$592.00.\$533.00 Current price is: \$533.00.



SKU: TAC-AT-KIT

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Tactee system is a set of assistive devices designed to restore autonomy and independence in the daily activities of people with mono or bilateral functional grasping deficits that aim to improve their quality of life. The Tactee system returns functional autonomy and independence in different daily activities such as eating, drinking, writing, etc.; eliminating the caregiver who's been considered indispensable so far.

When the grasping element is approached by a Tactee, tool, the magnetic force attracts it and allows a firm grip without using fingers. At the same time, its special design allows an easy release, permitting interchanging tools without any effort nor assistance, in total autonomy. The Tactee method involves the use of the system along the rehabilitation path of the users and promotes it in every context of everyday life, in an easy, decent and independent way. Tactee is a patented system based on an innovative magnetic coupling solution, between a grasping element worn by the hand's user and various interchangeable tools that guarantee high performance in autonomy and safety.

FEATURES

- **Ergonomic Design:** Handles are crafted to fit comfortably in the hand, reducing strain and enhancing grip, even for those with limited dexterity.
- **Modular Attachments:** A selection of interchangeable utensil heads, including spoons, forks, and knives, allows users to customize the cutlery to their specific needs.
- **Secure Connection:** The system features a simple yet secure mechanism for attaching and detaching the heads, ensuring they stay in place during use.



- **Durable Materials:** Constructed with high-quality, dishwasher-safe materials, the cutlery is easy to clean and designed to withstand daily use.
- **Non-Slip Handles:** The handles are textured or coated with a non-slip material to provide a stable grip, even for users with compromised hand strength.

TACTEE KIT CONTENTS

- **Grasping Element** the flexible strap made from polyurethane that makes it both hygienic and easy to close the grasping element using the mouth. To secure the grasping element further there is a thumb strap that prevents the adapter from sliding off the hand once loaded. Can be used on the right or left hand. Available in S, M, or L.
- Fork, Spoon & Knife Injection moulded and made of stainless steel and polycarbonate. The angle of the fork, spoon & knife can be adjusted by using a bending iron.
- Universal Adapter can be used with a toothbrush, razor, pen or other items.
- Bottle Holder has a velcro strap attached which makes it easy to attach a bottle or cup.
- **Pads** a set of dorsal pads that can be applied on the Grasping Element for users who are pressure-sensitive.
- Black Storage Case allows you to store, protect and transport the Tactee Kit contents safely.

SIZING

Available in 3 sizes. Measurement to be taken around MCP joints.

BROCHURE

View the Medimart Tactee brochure.