



TACTEE FORK

\$118.00



SKU: TAC-AT-TT02

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Tactee Fork is ideal for those with hand limitations. Use the Tactee Fork as an everyday aid for eating or cooking. Magnetic inserts allow users to quickly and easily connect and detach the fork to a [Tactee Grasping Element](#).

- Stable hold
- Injection molded and made of stainless steel and polycarbonate
- Angle of the fork can be adjusted by using a bending iron
- Dishwasher safe
- Suitable for right or left hand

Can also be used to practice functional gripping. Together with the Tactee Grasping Element, the fork supports you when eating and practicing the necessary movement sequences in functional gripping.

Accessories are sold separately or can be bought together in a [Tactee Kit](#) with a Grasping Element size of choice.