

TACTEE GRASPING ELEMENT

\$165.00



SKU: TAC-AT-TT01

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

Tactee is a new set of assistive devices with magnets that provides new opportunities to create a functional grip. Tactee consists of a Hand Grasping Element with three strong magnets that are placed next to each other with flexible joints between. The linear position of the magnets provides a secure and strong attachment. Designed to restore autonomy and independence in daily activities of people with mono or bilateral functional grasping deficits that aims to improve their quality of life. Suitable for neurological injuries, stroke, MS, Tetraplegia.

Tactee Grasping Element

The Tactee Grasping Element has a flexible strap made from polyurethane that makes it both hygienic and easy to close the grasping element using the mouth. To secure the grasping element further there is a thumb strap that prevents the adapter from sliding off the hand once loaded. Can be used on the right or left hand and paired with respective Tactee Accessories quickly and easily.

Sizing Guide:

Available in 3 sizes - small, medium and large. Measurement to be taken around MCP joints.

Use the Tactee Grasping Element with the following Tactee Tools & Accessories:

- [Tactee Fork](#)
- [Tactee Spoon](#)
- [Tactee Knife](#)
- [Tactee Universal Adapter](#)
- [Tactee Bottle Holder](#)
- [Tactee Pads](#)

- [Tactee Black Storage Case](#)

Accessories are sold separately or can be bought together in a [Tactee Kit](#) with a Grasping Element size of choice.