

TACTEE KNIFE

\$172.00



SKU: TAC-AT-TT04



PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Tactee Knife is ideal for those with hand limitations. Use the Tactee Knife as an everyday aid for eating or cooking. Magnetic inserts allow users to quickly and easily connect and detach the knife to a [Tactee Grasping Element](#).

- Stable hold
- Injection molded and made of stainless steel and polycarbonate
- Angle of the knife can be adjusted by using a bending iron
- Dishwasher safe
- Suitable for right or left hand

Accessories are sold separately or can be bought together in a [Tactee Kit](#) with a Grasping Element size of choice.