



TACTEE SPOON

\$118.00



SKU: TAC-AT-TT03

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The Tactee Spoon is ideal for those with hand limitations. Use the Tactee Spoon as an everyday aid for eating or cooking. Magnetic inserts allow users to quickly and easily connect and detach the spoon to a [Tactee Grasping Element](#).

- Stable hold
- Injection molded and made of stainless steel and polycarbonate
- Angle of the spoon can be adjusted by using a bending iron
- Dishwasher safe
- Suitable for right or left hand

Can also be used to practice functional gripping. Together with the Tactee Grasping Element, the spoon supports you when eating and practicing the necessary movement sequences in functional gripping.

Accessories are sold separately or can be bought together in a [Tactee Kit](#) with a Grasping Element size of choice.