



TECH-ASSIST ALL PURPOSE SLING

\$421.00 - \$433.00 Price range: \$421.00 through \$433.00



SKU: 23411

PRODUCT DESCRIPTION

PURPOSE & BENEFITS

The All Purpose Sling is a 2-in-1 sling designed to give greater flexibility of client use. It incorporates head support to allow for a full body lift but can also be used without the head support as a G.P. style sling by folding the head section out of the way. The head section foam can also be removed for comfort if leaving behind client. Webbing attachments on sling are covered to prevent chafing. Fully washable and quick drying. S.W.L 240KG.

SIZING